

DATA SHEET

Product Name: Tau-3G3 Antibody

Catalog #: TA-3012

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG2b_k

Supplied As: PBS pH 7.4 without additives

Dilution: 1 µg / ml Western Blot

Storage: at -20 °C (repeated thawing and freezing should be avoided)

Description: Tau 3G3 is a monoclonal antibody that reacts with Tau that is phosphorylated at amino acids 231 and 235 (T). Tau 3G3 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

For research use only. Not for use in humans.