

**DATA SHEET**

**Product Name:** Tau 8F10 Antibody

**Catalog #:** TA-3014

**Purity:** > 95% by SDS-PAGE

**Host/Isotype:** Mouse / IgG2b<sub>k</sub>

**Supplied As:** PBS pH 7.4 without additives

**Dilution:** 1 µg / ml Western Blot

**Storage:** at -20 °C (repeated thawing and freezing should be avoided)

**Description:** Tau 8F10 is a monoclonal antibody that reacts with all 6 isoforms of Tau at the Epitope: LADEVASLA. Tau 8F10 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

***For research use only. Not for use in humans.***