

DATA SHEET

Product Name: Tau E3 2B6 Antibody

Catalog #: TA-3018

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG2b_k

Supplied As: PBS pH 7.4 without additives

Dilution: 1 µg / ml Western Blot

Storage: At -20 °C (repeated thawing and freezing should be avoided)

Immunogen: Recombinant human Tau-441 in vitro aggregated

Specificity: Exon 3 of human Tau-441, amino acids 74-102

Reactivity: Tau-441, Tau-410

For research use only. Not for use in humans.